

Hazmat and Infection Control Catalog

**NIOSH Respirator
Standards Info**

N95 Respirators

**A, B and C Level
Protective Suits and Kits**

**Masks, Shields and
Protective Clothing**

**Gowns, Head Cover
Shoe and Boot Covers**

**Medical Gloves, ID Kits
and Lab Products**



**Recipient of the Vendor Excellence Award
from the U.S. Department of Defense**



Introduction and Catalog Index

Introduction

This catalog was developed to assist our existing customers as well as introduce our new customers to the products that many of you have requested over the years. Our goal is to develop closer ties and relationships with you by providing "One-Stop" shopping for your morgue, field and lab supply needs.

Although there are many common products offered herein, we have developed access to resources which allow us to locate and source almost any morgue or lab product commonly used in our field. Please call us or ask your Customer Service Representative for any products not included in this catalog offering. As recipients of the Department of Defense Vendor Excellence, we strive to provide only the most complete and responsive customer service in the industry.

This catalog as well as other information is available by contacting us toll free at (800) 345-2922. You may also download it from our internet site www.medicalproducts.com or email us at sales@medicalproducts.com.

Page Index

<u>Product Category</u>	<u>Page Number</u>
N95 and Cartridge Respirators	3 – 4
Tyvek ® Protective Suits and Coveralls	5 - 7
Tyvek ® Scene Responder Kits	7 - 9
Procedure Masks and Shields	9 - 10
Head Cover, Gowns, Shoe Covers, Boot Covers, Gowns, Lab Coats, Sleeve Protectors and Scrubs	10 - 14
Medical Gloves	14 – 16
Lab Supplies, ID Products and Kits	17 – 18
NIOSH Standards, Sizing Charts and Mat'l Information	19 – BACK CVR



NIOSH and CDC Approved N95 Respirator

- Magic Arch® polypropylene strips encased in mask material to provide a comfortable breathing chamber.
- Aluminum strips on both top and bottom to ensure a custom fit and seal.
- Fluid resistant outer facing with cool, soft inner material.
- Mask material is fiberglass and latex free with comfortable dual knit elastic headbands.
- NIOSH Approved Universal "One-Size-Fits-All" Design. No need to stock multiple size masks for each user. NIOSH Cert: TC-84A-0427
- High filtration efficiency with a BFE at over 99% at .1 micron
- Packed 35 Units Per Box/ 6 Boxes Per Case



Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-N95-695	UNIV	TEAL/WHT	210	6	2.72
MASK-FIT-TEST	QUAL FIT TEST KIT		1	1	.45



Gerson N95 Cone Style Respirator

- NIOSH N95 Approved Particulate Respirator
- Lightweight, comfortable and easy to use, staple free construction
- Meets CDC guidelines for TB exposure control
- Fluid resistant
- Durable double-shell resists crushing & provides durability in hot, humid conditions



Part No	Size	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-N95-1730	UNIV	Cone Style N95 Respirator	WHT	240	6	2.72



Disposable "Odor Eater" OV Cone Respirator

- Conforms to NIOSH 42 CFR Part 84, N95 class
- Active Carbon Melt Blown (ACMB) fibers are resin free and provide excellent protection against organic vapors and welding fumes
- Filtration 95% effective against particulates and aerosols free of oil
- Adjustable plastic nose piece and inner foam provide a comfortable facial seal



Part No	Description	Size	Color	Qty/Cs	Lb/Cs	Kg/Cs
MASK-N95-87120V	OV Cone Style N95 Mask w/Ex Valve	UNIV	WHT	144	7	3.18

We encourage you to check inside the rear cover for NIOSH Specifications

Ask your Customer Service Representative about our Fit Test Kits



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200

www.medicalproducts.com • sales@medicalproducts.com

Gerson Reusable Cartridge Respirator



- Low maintenance and low cost cartridge respirator
- NIOSH Approved
- Soft, non-allergenic face piece
- Replaceable cartridges and filters sold separately

Part No	Size	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-GER-9100	S/M/L	HALF MASK RESP	BLK	24	6	2.72
P100, Organic Vapor and/or Acid Gas Filter Cartridges Available						

North Safety Half Face Cartridge Respirator



- Silicone reusable NIOSH half face respirator protection
- Lightweight, comfortable and easy to use
- Available in three (3) sizes for custom facial fit
- Filter cartridges available for most particle and vapor risks
- Silicone material to assure long term use and a positive facial seal (Also available in less expensive elastomeric material)

Part No	Size	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-NOS-7700	S/M/L	HALF MASK RESP	BLK	1	2	.9
P100, Organic Vapor and/or Acid Gas Filter Cartridges Available						

North Safety Full Face Cartridge Respirator



- Reusable NIOSH full face respirator protection
- Lightweight, comfortable and easy to use
- Available in three (3) sizes for custom facial fit
- Filter cartridges available for most particle and vapor risks
- Silicone material to assure long term use and a positive facial seal

Part No	Size	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-NOS-7600	S/M/L	FULL FACE RESP	BLK	1	4	1.81
P100, Organic Vapor and/or Acid Gas Filter Cartridges Available						

We encourage you to check inside the rear cover for NIOSH Specifications




Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200

www.medicalproducts.com • sales@medicalproducts.com

Level "C" Protective Suits



When you see this new label, you'll know it's New Tyvek®!



Tyvek ® Basic Coverall

- Tyvek ® Coverall w/ Elastic Wrist and Ankles or set in sleeves
- Sewn seam construction with full front nylon zipper
- Semi-impervious "New" Tyvek ® Material
- Offers a dry particulate barrier and some light splash protection

Set In Sleeve Basic Coverall (Meets ANSI/ISEA 101-1996 specifications)

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CVRALL-1412-XX	SEWN	WHITE	25	8	3.63

Same as above but with Elastic Wrist and Ankles

CVRALL-1417-XX	SEWN	WHITE	25	8	3.63
----------------	------	-------	----	---	------

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X **



Tyvek ® Coverall With Hood and Boots

- One piece Tyvek ® Suit w/ Attached Hood and Anti-Skid Boots
- Sewn seam construction with full front nylon zipper
- Semi-impervious, offers a dry particulate barrier and some light splash protection

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CVRALL-1414-XX	SEWN	WHITE	25	10	4.54

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X **

Level "B" and "C" Protective Suits

Tychem QC ® Coverall With Hood and Boots



- One piece Tychem ® QC Suit w/ Attached Hood and Boots
- Available with sewn or sealed seams
- Full front nylon zipper, elastic wrists and face
- Impervious to blood borne fluids as well as many acids and alkalis
- Can provide protection from blood borne pathogens, Level "B" Hazmat Protection from several chemical agents in dry aerosol form, light splash protection for sewn version and heavy splash protection in the sealed seam construction.

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CVRALL-5414-XX	SEWN	YELLOW	12	7	2.20
CVRALL-70150-XX	SEALED	YELLOW	6	8	3.63

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X **

**** SEE BACK COVER OF THIS CATALOG FOR SIZING CHART AND CONSTRUCTION SPECS**



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200

www.medicalproducts.com • sales@medicalproducts.com



Tychem SL ® Coverall With Hood and Boots

- One piece Tychem SL ® Suit w/ Attached Hood and Boots w/Straps
- Full front "storm flapped" zipper, with elastic wrists and face openings
- Impervious Tyvek ® spun bonded olefin by Dupont laminated with Saranex ® 23-P film from Dow Chemical
- Available with bound or sealed seams
- Can provide Level "B" Hazmat Protection in the presence of most hazardous chemicals and bio-pathogens.

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CVRALL-44414-XX	BOUND	WHITE	12	15	6.8
CVRALL-72150-XX	SEALED	WHITE	6	8	3.63

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X **



Tychem BR ® Coverall With Hood and Boots

- One piece Tychem BR ® Suit w/ Attached Hood and Boots
- Full front "storm flapped" zipper, with elastic wrists and face openings
- Newest Dupont high level material for Evaluating Protective Clothing Materials (ASTM F1001) for over 200 chemicals and agents
- Sealed seam construction
- Provides solid Level "B" Protection for most Hazmat exposures

Part No.	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CVRALL-BR150-XX	SEALED	YELLOW	6	8	3.7

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X **



Tychem BR ® Level "B" Encapsulated Suit

- One-piece Tychem BR ® Encapsulated rear entry suit
- 48" zipper, double storm flap with Velcro®, 20 mil PVC face shield, elastic wrists, 2 exhaust ports with shrouds, attached sock boots with boot flaps and expanded back
- Newest Dupont high level material for Evaluating Protective Clothing Mat'ls (ASTM F1001) for over 200 chemicals and agents.
- Sealed seam construction including face shield
- Provides Level "B" Protection for most HAZMAT exposures

Part No.	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CVRALL-BR450-XX	SEALED	LT GRN	3	13	5.9

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X **

**** SEE BACK COVER OF THIS CATALOG FOR SIZING CHART AND CONSTRUCTION SPECS**



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



Tychem TK ® Level "A" Encapsulated Suits

- Fully encapsulated front or rear entry vapor proof suit with expanded back to accommodate SCBA gear
- Sealed seams inside and out with 48" covered zipper
- Double storm flap with Velcro ® closure
- Two layer (10 mil Teflon/40 mil PVC) face-shield
- Three layer knit glove combination of Butyl, North Silvershield and outer Kevlar material
- Dupont's highest level material for protective suits

Part No	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
BIOSUIT-TK640	Front Entry Tychem Level "A" Suit	LT GRN	1	10	3.18
BIOSUIT-TK650	Rear Entry Tychem Level "A" Suit	LT GRN	1	10	3.18

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X/5X **

Level "B and C" Protective Suit Kits



New Tyvek ® Disaster and Crime Scene Response Kit

- Tyvek ® Coverall w/Hood, Boots, Elastic Wrist and Face
- 2 Pair of Large Nitrile Gloves
- Pair of Safety Goggles
- Packaged in Ziploc reusable poly bag
- N95 NIOSH Respirator
- Bio-Hazard Disposal Bag

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
KIT-IFC-01-XX	SEWN	WHITE	10	9	4.4

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X/5X **



Tychem QC ® Disaster and Crime Scene Response Kit

- Tychem QC ® Coverall w/Hood, Boots, Elastic Wrist and Face
- 2 Pair of Large Nitrile Gloves
- Pair of Safety Goggles
- Packaged in Ziploc reusable poly bag
- N95 NIOSH Respirator
- Bio-Hazard Disposal Bag

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
KIT-IFC-01QC-XX	SEWN	YELLOW	10	15	6.8

XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X/5X **



All Protective Suit Kits Available With Rugged Gear Carry and Storage Bag !



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



Tychem QC ® Disaster and Crime Scene Response Kit

- Tychem QC ® Coverall w/Hood, Boots, Elastic Wrist and Face
- 2 Pair of Large Nitrile Gloves
- Pair of Safety Goggles
- Packaged in Ziploc reusable poly bag
- N95 NIOSH Respirator
- Bio-Hazard Disposal Bag

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
KIT-IFC-01QC-XX	SEWN	YELLOW	10	15	6.8
XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X/5X **					



Tychem SL ® Disaster and Crime Scene Response Kit

- Tychem SL ® Coverall w/Hood, Boots, Elastic Wrist and Face
- 2 Pair of Large Nitrile Gloves
- Pair of Safety Goggles
- Packaged in Ziploc reusable poly bag
- N95 NIOSH Respirator
- Bio-Hazard Disposal Bag

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
KIT-IFC-01SL-XX	SEALED	WHITE	3	8	3.63
XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X/5X **					



Tychem SL ® HD Chemical Agent Disaster Response Kit

- Tychem SL ® Coverall w/Hood, Boots, Elastic Wrist and Face
- 2 Pair of Large Nitrile Gloves
- 10 yd. Roll of Grey Tape
- Pair of 12" High Hazmat Ribbed Sole Yellow Boot Covers
- Packaged in Ziploc reusable poly bag
- Pair of Chemical Res. Gloves
- Bio-Hazard Disposal Bag

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
KIT-HZMT-72150-XX	SEALED	WHITE	3	8	3.63
XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X/5X **					



Tychem BR ® HD Chemical Agent Disaster Response Kit

- Tychem BR ® Coverall w/Hood, Boots, Elastic Wrist and Face
- 2 Pair of Large Nitrile Gloves
- 10 yd. Roll of Grey Tape
- Pair of 12" High Hazmat Ribbed Sole Yellow Boot Covers
- Packaged in Ziploc reusable poly bag
- Pair of Chemical Res. Gloves
- Bio-Hazard Disposal Bag

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
KIT-HZMT-BR150-XX	SEALED	WHITE	3	8	3.63
XX – Available Sizes: SM/MED/LG/XL/2X/3X/4X/5X **					



All Protective Suit Kits Available With Rugged Gear Carry and Storage Bag !



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200

www.medicalproducts.com • sales@medicalproducts.com



Post Decontamination Cover Up Kit

1. Little White Reusable Plastic Ziploc Bag
2. Bio-Hazard Disposal Bag
3. Custody Card
4. Bio-Hazard Disposal Bag

Kit No. _____
Kit Date _____

Kit Issued To: _____
Kit Issued By: _____
Kit Issued At: _____

Kit Returned To: _____
Kit Returned By: _____
Kit Returned At: _____

Kit Returned Date: _____

Kit Returned Status: _____



Post Decontamination Cover-up Suit Kit

- Comfortable light blue SMS coverall with attached hood and boots
- Bio-Hazard Disposal Bag
- Custody card and polybag for storage and transport of personal items such as jewelry, watches, wallets, etc.
- Packaged in Ziploc reusable poly bag

Part No	Seam Const.	Color	Qty/Cs	Lbs/Cs	Kg/Cs
KIT-HZMT-01*	SEWN	BLUE	10	7	3.18
XX – Available Sizes: SM/MED/LG/XL/2X/3X**					



Isolation Ear Loop Procedure Mask

- Polypropylene .3 micron isolation mask
- Comfortable ear loop procedure mask design
- Full width aluminum nose bridge band for good facial seal

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-BL5005	UNIV	BLU	500	5	2.27



Surgeon's Face Mask w/Ties

- Polypropylene .3 micron isolation mask
- Tie back mask design w/ 7" ties
- Full width aluminum nose bridge band for good facial seal

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-BL5001	UNIV	BLU	300	6	2.72



Isolation Mask and Shield Combination

- Polypropylene 3 micron isolation mask with attached 7 mil Anti-fog eye shield
- Comfortable ear loop procedure or tie back mask design
- Full width aluminum nose bridge band for good facial seal
- Full "wrap around" shield design for maximum eye protection

Part No	Style	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-SM2005X	EARLOOP	BLUE	100	5	2.27
MASK-BL8101S	TIES	BLUE	100	5	2.27



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200

www.medicalproducts.com • sales@medicalproducts.com



Cone Style Gen Purpose Mask

- Cone style mask w/elastic head band
- Aluminum nose bridge clip for good facial seal
- Better than 95% BFE filtration protection
- Easy to don and comfortable to wear

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-BL2020	UNIV	BLU	300	5	2.28



Anti-Fog Face Shields

- Premium 7 mil optical grade anti-fog film face shield material
- Soft foam forehead perspiration cushion
- Soft elastic comfort band for comfortable long-term use or adjustable Velcro style head band
- Conforms to OSHA and CDC recommended facial protection

Part No	Headband Style	Color	Qty/Cs	Lbs/Cs	Kg/Cs
MASK-2807	COMFORT BAND	CLR	100	10	4.36
MASK-2802	VELCRO BAND	CLR	100	10	4.36



Spunbond Bouffant Caps

- Spunbond poly-pro material w/non-binding elastic headband
- Comfortable and cool to wear
- Provides splash or mist fluid protection

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CAP-120070-SB-21	21" (53 cm)	BLU	1,000	16	7.26
CAP-120050-SB-21	21" (53 cm)	WHT	1,000	16	7.26
CAP-120070-SB-24	24" (61 cm)	BLU	1,000	16	7.29
CAP-120050-SB-24	24" (61 cm)	WHT	1,000	16	7.29



Tyvek ® Bouffant Cap

- Tyvek ® material w/non-binding elastic headband
- Provides greater fluid protection than the spunbond bouffant caps above

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CAP-80221-TVK	21" (53 cm)	WHT	500	7	3.2
CAP-80224-TVK	24" (61 cm)	WHT	500	7	3.2



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



Surgeon Cap

- Scrim reinforced fabric material w/ties
- Comfortable splash and mist protection

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
CAP-SURG-B	UNIV	BLU	500	13	5.9



Poly-Film Thumb Hole Gown

- All poly film impervious material
- Unique thumb loop sleeve extension design eliminates sleeve riding up from under exam glove and exposing bare skin to pathogen environment
- Easy "slip-on" boat-neck style collar for easy donning and open back design for and comfort during long procedures



Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
GOWN-8572	UNIV	BLU	100	10	4.6



Poly Isolation Gown w/Elastic Knit Cuffs

- Open back impervious isolation gown with cotton "Comfort Cuff"
- Elastic sleeves and below knee length

Part No	Color	Qty/Cs	Lbs/Cs	Kg/Cs
GOWN-110070-XX	BLUE	50	9	3.8

XX – Specify Size – LG, XL or 2X



Poly Lab Apron

- 1.75 mil polyethylene general purpose lab apron
- Apron has neck strap and waist ties

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
APRON-52846-175	28 X 46"	WHITE	500	30	13.6



SMS Comfort Cool Lab Jacket

- Comfortable elastic knit cotton collar and cuffs
- Secure plastic snap front closure with dual pockets
- Tri-Layer Melt Blown Spun bond material allowing for cooler wear yet still having excellent impervious properties

Part No	Color	Qty/Cs	Lbs/Cs	Kg/Cs
LJ-SMS-XX	BLUE	30	11	4.99

XX – Specify Size - SM, MED, LG, XL or 2X



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



SMS Comfort Cool Lab Coat

- Traditional collar with comfortable elastic cotton knit cuffs
- Secure plastic snap front closure with dual pockets
- Tri-Layer Melt Blown Spun bond material allowing for cooler wear yet still having excellent impervious properties

Part No	Color	Qty/Cs	Lbs/Cs	Kg/Cs
LC-SMS-MED	BLUE	50	12	5.44
XX – Specify Size - SM, MED, LG, XL or 2X				



Spun bond Disposable Scrubs

- Comfortable fitting semi-impervious spun bond material
- Pants have draw string waistband with pocket

Part No	Color	Qty/Cs	Lbs/Cs	Kg/Cs
SUIT-11005-XX	BLUE	50	27	12.25
XX – Specify SM, MED, LG, XL or @X				



Tyvek® Aprons and Sleeve Protectors

- Apron is 28" Wide x 36" Long with sewn ties
- Aprons packed 100 units per case
- Sleeves are 18" long and have elastic top, cuffs and sewn seams.
- Sleeves Packed 200 Units (100 Pair) Per Case

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
APRON-601	28 X 36"	WHITE	100	6	2.91
SLEV-85018	18"	WHITE	200	7	3.4



Tyvek® QC Aprons and Sleeve Protectors

- Apron is 28" Wide x 36" Long with sewn ties
- Aprons packed 100 units per case
- Sleeves are 18" long and have elastic top, cuffs and sewn seams.
- Sleeves Packed 200 Units (100 Pair) Per Case

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
APRON-650	28 X 36"	YELLOW	100	9	4.4
SLEV-86018	18"	YELLOW	200	7	3.4



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



Anti-Skid Shoe Covers

- Spun bond Polypropylene Semi Impervious Material
- Anti-Skid rubber sole for good traction and workplace safety
- Univ size packed 300 units (150 Pair) per Case
X-Large size packed 200 units (100 pair) per case

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
SHOECVR-110070-UNI	UNIV	BLUE	300	6	2.92
SHOECVR-110070-XL	X-LARGE	BLUE	200	7	3.18

Tyvek ® Shoe Covers

- Elastic top shoe cover made of Dupont Tyvek ®
- Semi-impervious for wet or dry floor use
- Packed 400 Units (200 Pair) Per Case



Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
SHOECVR-901-UNIV	SM/MED	WHITE	400	7	3.18
SHOECVR-901-XL	LG/XL	WHITE	400	8	3.63

Tyvek ® Anti-Skid Shoe Covers

- Elastic top shoe cover made of Dupont Tyvek ® Anti-Skid Mat'l
- Semi-impervious for wet or dry floor use
- Packed 400 Units (200 Pair) Per Case

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
SHOECVR-901NS-UNIV	SM/MED	GREY	400	7	3.18
SHOECVR-901NS-XL	LG/XL	GREY	400	8	3.63

Tyvek ® Boot Covers

- 15" Tyvek ® Elastic top boot cover made of Dupont Tyvek ®
- Semi-impervious for wet or dry floor use
- Packed 400 Units (200 Pair) Per Case



Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
BOOTCVR-903-SMED	SM/MED	WHITE	400	10	4.54
BOOTCVR-903-LGXL	LG/XL	WHITE	400	10	4.54

Tyvek ® Anti-Skid Boot Covers

- 15" High elastic top made of Dupont Tyvek ® Anti-Skid Material
- Semi-impervious for wet or dry floor use
- Packed 400 Units (200 Pair) Per Case

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
BOOTCVR-903NS-SMED	SM/MED	GREY	400	10	4.54
BOOTCVR-903NS-LGXL	LG/XL	GREY	400	10	4.54



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



Aqua Trak® Impervious Shoe Cover

- High quality poly-coated spun bond material
- Seams are heat sealed to make this shoe cover as impervious as possible while still offering safe, comfortable protection.
- Packed 150 Units (75 Pair) Per Case

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
SHOECVR-AT11111-B	UNIV	BLUE	150	6	2.92
SHOECVR-AT21112-B	X-LARGE	BLUE	150	7	3.18



Aqua Trak® Impervious Boot Cover

- 18" High quality poly-coated spun bond material
- Seams are heat sealed to make this shoe cover as impervious as possible while still offering safe, comfortable skid resistant workplace protection.
- Packed 200 Units (100 Pair) Per Case

Part No	Size	Color	Qty/Cs	Lbs/Cs	Kg/Cs
BOOTCVR-BT91012B	UNIV	BLUE	200	14	6.35
BOOTCVR-BT91013B	X-LARGE	BLUE	200	14	6.35



HazMat Boot Covers

- 12" High seamless latex boot covers made of .70mm latex rubber for extra protection in hostile chemical environments
- "Hi visibility" yellow color with ribbed sole for non-slip protection

Part No	Size (Shoe)	Color	Qty/Cs	Lbs/Cs	Kg/Cs
BOOTCVR-97591-L	LG (9-11)	YELLOW	50 (25 PR)	24	10.9
BOOTCVR-97591-XL	XL (10-15)	YELLOW	50 (25 PR)	24	10.9



Silvershield® Chemical Resistant Glove

- Made of Norfoil®, a lightweight, flexible laminate material
- Highly impermeable to more chemicals and solvents than any other type of glove on the market today
- Unique, thin and flexible film that conforms to the hand for greater dexterity
- Ambidextrous (Avoids Hand Pairing)
- Can be used as an inner glove thus allowing maximum protection in heavy-duty applications where dangers of mechanical damage are high.

Part No	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-NOR-SSG- **	Silvershield® Chemical Resistant Gloves	GRY	100	4	1.8

** Available Sizes – 6, 7, 8, 9, 10 and 11



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



Latex Exam Gloves

- 7 Mil Palm – 5 mil Cuff White Latex Exam Glove
- Available both lightly powdered and powder free
- Packed 100 per box - 10 Boxes Per Case

Part No	Style	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-DC-LP- *	9" EXAM GLV LT POWDER	1,000	20	9.07
GLV-DC-PF- *	9" EXAM GLV PDR FREE	1,000	20	9.07
GLV-OGP-PF- *	9" EXAM GLV PF TEX PALM	1,000	20	9.07

* SPECIFY SIZE – SM, MED, LG OR XL



BarrierMax™ Long Cuff Latex Exam Glove

- White High-Risk Long Cuff Powder Free Latex Exam Glove
- 12 Mil Palm – 8 mil 12" cuff w/textured grip-tight fingers
- Packed 50 per box – 10 Boxes Per Case

Part No	Style	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-BMAX-PF *	BARRIERMAX PF GLV	500	20	9.07

* SPECIFY SIZE – SM, MED, LG OR XL



ICU911™ EMS High Risk Long Cuff Glove

- High-Risk Long Cuff Powder Free Latex Exam Glove
- 10 Mil Palm – 6 mil 12" cuff with textured grip-tight fingers
- Packed 50 Per Box – 10 Boxes Per Case

Part No	Style	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-ICU911-PF-*	ICU911 PF LATEX GLV	500	20	9.07

* SPECIFY SIZE – SM, MED, LG OR XL



HaloKote HD Chemotherapy Glove

- Green High-Risk Toxic Chemical Resistant Glove
- 15 mil textured fingers / 13 Mil Palm / 8 Mil 12" Cuff
- Packed 50 per box – 10 Boxes Per Case

Part No	Style	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-DSK-PF- *	DSK POWDER FREE	500	22	10.0

* SPECIFY SIZE – SM, MED, LG OR XL

GLOVE SIZING CHART

GLOVE SIZE (US)	6-7	7-8	8-9	9 AND UP
EXAM GLV SIZE	SMALL	MED	LARGE	EXTRA LARGE



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



Frontline™ Nitrile Glove

- Blue High-Risk Powder Free Nitrile Exam Glove
- 6/5 mil Palm 9" 4 mil cuff with textured grip-tight fingers
- Packed 100 per box – 10 Boxes Per Case

Part No	Style	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-FL9-PF *	9" CUFF NITRILE EXAM GLV	1,000	20	9.07
* SPECIFY SIZE – SM, MED, LG OR XL				



Frontline™ Long Cuff Nitrile Glove

- Blue High-Risk Long Cuff Powder Free Nitrile Exam Glove
- 6/5 Mil Palm 12" 4 mil cuff with textured grip-tight fingers
- Packed 100 per box – 10 Boxes Per Case

Part No	Style	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-FL12-PF *	12" CUFF NITRILE EXAM GLV	1,000	22	9.98
* SPECIFY SIZE – SM, MED, LG, XL or 2XL (2XL 900/CS)				



Defender 7 Mil Powder Free Nitrile Glove

- Slate Color High-Risk Powder Free Nitrile Exam Glove
- 8 mil textured fingers / 7 Mil Palm / 5 Mil Long Cuff

Part No	Style	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-DFR10-*	10" CUFFNITRILE GLOVE	1,000	23	10.43
GLV-DFR12-*	12" CUFF NITRILE GLOVE	500	16	7.3
* SPECIFY SIZE – SM, MED, LG OR XL				



Safeknit Spectra® Cut Resistant Glove

- Soft and Comfortable Spectra® Fiber Material
- Excellent cut resistant properties yet allows good tactile feel
- Stands up to laundering even with bleach !
- Ambidextrous - worn as a liner under exam or autopsy gloves



Part No	Style/Color	Qty/Cs	Lbs/Cs	Kg/Cs
GLV-75025- *	KNIT CUT RESIST WHITE	12	1	.45
* SPECIFY SIZE – 6, 7, 8, 9, 10, 11				

GLOVE SIZING CHART				
GLOVE SIZE (US)	6-7	7-8	8-9	9 AND UP
EXAM GLV SIZE	SMALL	MED	LARGE	EXTRA LARGE



Medicalproducts Ltd., Inc.


Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com



Spinal and Body Fluid Recovery Syringe Needles

- 304 Surgical Stainless Steel w/Standard Luer-Lok Hub
- 13 GA Long Length Fluid Recovery Size w/metal hub available

Part No	Size	Color	Qty/Bx	Lbs/Bx	Kg/Bx
PM-BD405181	22 GA X 3.5" (8.9 CM)	SS	25	.5	.23
PM-BD305196	18 GA X 1.5" (3.8 CM)	SS	100	.5	.23
PM-BD405184	18 GA X 3.5" (8.9 CM)	SS	25	.5	.23
PM-BD305198	16 GA X 1.5" (3.8 CM)	SS	10	.5	.23
PM1622-1013R	12 GA X 4.5" (11.5 CM)	SS	12	.5	.23



Disposable Syringes With Luer Tips

- High quality plastic syringes with Luer Tips
- Individual sterile packaging

Part No	Capacity	Color	Qty/Bx	Lbs/Bx	Kg/Bx
SYRINGE-510125	10 CC	CLR	100	1	.45
SYRINGE-520125	20 CC	CLR	25	.5	.23
SYRINGE-530125	30 CC	CLR	25	.6	.28
SYRINGE-560125	60 CC	CLR	25	.7	.32



Scalpel Blades

- Surgical quality scalpel blades
- Pkg. sterile 100 blades per box

Part No	Size	Material	Qty/Bx.	Lbs/Bx.	Kg/Bx.
SBLADE-22-C	22 GA	CARBON STEEL	100	.5	.23
SBLADE-22-SS	22 GA	STAINLESS	100	.5	.23
SBLADE-23-C	23 GA	CARBON STEEL	100	.5	.23
SBLADE-23-SS	23 GA	STAINLESS	100	.5	.23
SBLADE-60-C	60 GA	CARBON STEEL	100	.5	.23
SBLADE-60-SS	60 GA	STAINLESS	100	.5	.23
SBLADE-70-C	70 GA	CARBON STEEL	100	.5	.23

Scalpel Blade Handles

- Surgical quality stainless steel material
- Available in both solid and hollow handle design



Solid Handle Design

Part No	Size	Blade Sizes	Qty/Bx	Lbs/Bx	Kg/Bx
SBHAND-04-SS	4" (10.2 CM)	20,21,22,23,24,25	1	.5	.23
SBHAND-04L-SS	7" (17.8 CM)	20,21,22,23,24,25	1	.5	.23

Hollow Handle Design


Part No	Size	Blade Sizes	Qty/Cs	Lbs/Cs	Kg/Cs
SBHAND-06-SS	6" (15 CM)	20,21,22,23,24,25	1	.5	.23
SBHAND-08-SS	6" (15 CM)	60,70	1	.5	.23



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200

www.medicalproducts.com • sales@medicalproducts.com



"Easy Write On" Metal ID Tag

- Easy Write-On metal tag with attachment wires
- Ideal for use in extreme environmental conditions
- Soft metal allows use of almost any marking or writing implement including sticks, keys etc.

Part No	Packaging	Qty/Cs	Lbs/Cs	Kg/Cs
TAG-METAL-01	METAL TAG ONLY	100	1	.45

+

Serial Numbered Tamper Resistant Body Bag Seal

- 8" Long serial numbered plastic ID tag
- Sequentially numbered in lots of 100
- Designed to fit Medicalproducts body bag zippers
- Custom number series available (1,000 pcs min.)



Part No	Qty/Cs	Lbs/Cs	Kg/Cs
TAG-MPSEAL	100	1	.45



Write-On Wrist or Ankle ID Band

- Write-On band with "Use Once" permanent snap style closure
- Each case packed with marking pen
- Available in vinyl or polypro material

Part No	Material	Qty/Cs	Lbs/Cs	Kg/Cs
MC-01-4440-W	WHITE VINYL	250	2	.90
MC-01-5040-W	WHITE POLY-PRO	500	3	1.36



Body Bag Identification Kit

- Re-sealable, self adhesive Ziploc pouch for paperwork
- Red Bio-Hazard Tag
- Yellow Contamination Tag
- 3 White tags for toe, exterior of pouch & personal effects ID

Part No	Qty/Cs	Lbs/Cs	Kg/Cs
BBID-KIT	50	1	.45



Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com

NIOSH Respirator Filter Standards



1. Respirators with N100, N99, and N95 filters (99.97%, 99%, and 95% efficient filters) may be used for any solid or non-oil containing particulate contaminant.
2. Respirators with R100, R99, and R95 filters (99.97%, 99%, and 95% efficient filters) may be used for any particulate contaminant. If used for an oil containing particulate, a one- shift (8 hour) use limit applies.
3. Respirators with P100, P99, and P95 filters (99.97%, 99%, and 95% efficient filters) may be used for any particulate contaminant.

All NIOSH certified respirators have an assigned "protection factor". The "protection factor" indicates how much the respirator will reduce a wearer's exposure to airborne particles assuming a good fit.

The goal of using a respirator is to reduce the exposure to an acceptable level that will not adversely affect the wearer. The assigned protection factor predicts the potential reduction in exposure that would occur if one were properly wearing that type of respirator properly fit to the face.

Protection Level Table

Concentration of Airborne Particles	Respirator Protection Factor	Protection Level
10,000 Part. Per Cubic Foot	10	This means the person could expect to reduce their exposure by 10 times resulting in a possible inhalation of 1000 particles.
	50	This means the person could expect to reduce their exposure by 50 times resulting in a possible inhalation of 200 particles.
	1000	This means the person could expect to reduce their exposure by 1000 times resulting in a possible inhalation of 10 particles.
	10,000	This means the person could expect to reduce their exposure by 10,000 times resulting in a possible inhalation of 0 particles.


 Level A	Level A represents the greatest danger of respiratory, eye or skin damage from hazardous vapors, gases, particulates, sudden splash, immersion or contact with hazardous materials. It calls for total encapsulation in a vapor tight chemical suit with self-contained breathing apparatus (SCBA) or supplied air and appropriate accessories.
 Level B	Level B situations call for the highest degree of respiratory protection but a lesser need for skin protection. It calls for SCBA or positive pressure supplied air respirator with escape SCBA, plus hooded chemical resistant clothing (overalls and long sleeved jacket; coveralls; one or two piece chemical-splash suit; or disposable chemical-resistant coveralls boots and gloves.
LEVEL C	Level C is the lowest level of hazmat protection. Level "C" requires NIOSH N95 or higher level of respiratory protection but a minimal need for skin protection. It calls for half mask or full-face negative pressure respirators or powered air respirators (PAPR's). Level "C" requires use of a hooded chemical resistant suit or two piece chemical resistant suit, chemical resistant boots and gloves. Respirator must be able to remove <u>all</u> chemical and biological agents and the agent must pose <u>no splash hazard</u> .




Medicalproducts Ltd., Inc.

Toll Free (800) 345-2922 • Fax (800) 372-5649 • Ph: (678) 413-9200
www.medicalproducts.com • sales@medicalproducts.com


ANSI Coverall Sizing Chart



Encapsulated Suit Sizing Chart




Seam Construction for Protective Clothing




Serged Seam

Exposure limited to dry particulate agent or light splash of liquids not harmful to skin.



Sewn and Bound Seam

Higher degree of risk of exposure to agents in either liquid or dry form.



Heat Sealed Seam

Highest Degree of protection from harmful liquid, dry chemical or unknown agents

DISCLAIMER

THE INFORMATION CONTAINED IN THIS CATALOG IS PROVIDED TO YOU "AS IS," FOR YOUR INTERNAL INFORMATIONAL PURPOSES ONLY, WITHOUT ANY REPRESENTATION OR WARRANTY OF ACCURACY OR COMPLETENESS OF INFORMATION OR OTHER WARRANTY OF ANY KIND, INCLUDING ANY IMPLIED WARRANTY OF QUALITY, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT.

IN NO EVENT WILL MEDICALPRODUCTS LTD., INC., ITS AGENTS OR REPRESENTATIVES BE LIABLE TO ANY PARTY FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES FOR USE OF THIS CATALOG OR RELIANCE UPON ANY INFORMATION OR MATERIAL ACCESSED VIA IT OR ANY OTHER REFERENCED SOURCE INCLUDING, BUT NOT LIMITED TO, DAMAGES ARISING FROM LOSS OF PROFITS, BUSINESS INTERRUPTION, LOSS OF DATA OR PERSONAL INJURY EVEN IF MEDICALPRODUCTS LTD., INC., ITS AGENTS OR REPRESENTATIVES ARE EXPRESSLY ADVISED ABOUT THE POSSIBILITY OF SUCH DAMAGES, TO THE FULLEST EXTENT ALLOWABLE BY LAW. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OF IMPLIED WARRANTIES AND REPRESENTATIONS, SO THE ABOVE EXCLUSIONS MAY NOT APPLY TO YOU.



Medicalproducts Ltd., Inc.

Address: 341 Bullis Street, Eagle Pass TX 78852
 Mail Address: P.O. Box 80685, Conyers GA 30013-0685
 Toll Free Sales (800) 345-2922 • Fax: (800) 372-5649
 Office: (678) 413-9200 • Fax: (770) 860-8845

Internet: www.medicalproducts.com Email: sales@medicalproducts.com

Copyright © 2002, 2009, 2010 • Medicalproducts Ltd., Inc. • All Rights Reserved