

FirstCall™ Mortuary Products



**Recipient of the Department of
Defense Award for
Vendor Excellence**

Introducing our *NEW Line of FirstCall™ Mortuary Products*



All FirstCall™ Products
Made In U.S.A.



Magna Ultra Tuff © FirstCall™ Pouch

Full length pouch design with openings for webbing attachment to the cot and HD Nylon Zipper w/dual pulls

Made of Magna Ultra Tuff © Heavy Duty Stain Resistant Nylon Fabric coated with polyurethane for moisture protection

420 Denier nylon inner liner to further enhance durability

Document pouch with Velcro closure to protect paperwork

Designed to fit standard mortuary cots like our FCALL-COT-100 (as shown)

Available in Burgundy (-85), Navy Blue (-71) (shown), Hunter Green (-78) and Black (-60)

Part No	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
FCALL-200-XX	24 x 74 x 14" Ultra Tuff FirstCall™ Pouch	XX	1	6	2.3

FirstCall™ Cover

20 x 76 x 14 available in our brushed plush poly fabric or our Magna Ultra Tuff © Cordura

Polyurethane backing for easier cleaning

With or without elastic at the ends

Customized embroidery on one or both sides

Document pouch with Velcro closure

Available in Burgundy (-85), Navy Blue (-71) and Hunter Green (-78)



Part No	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
FCALL-CVR-XX	20 x 76 x 14" Brushed Poly FirstCall™ Cover	XX	1	2	.8



FirstCall™ Church Truck Drape

23 x 64 x 23 Brushed Fleece Nylon Fabric

Fits most church trucks

Even pleating with reinforced holes

Available in FirstCall™ Colors listed above

Part No	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
FCALL-DRAPE-XX	23 x 64 x 23" Brushed Poly FirstCall™ Drape	Burg	1	3	1.5



Medicalproducts Ltd. Eagle Pass TX 78852

Toll Free: (800) 345-2922 Fax: (800) 372-5649 Ph: (678) 413-9200

Internet Site: www.medicalproducts.com Email: sales@medicalproducts.com

A Service Disabled Veteran Owned Company

Introducing our *NEW FirstCall™* Multi-Level One Man Mortuary Cot



**All FirstCall™ Products
Made In U.S.A.**

- Multi-level locking positions and roll-in style allows one person to easily perform retrieval
- Special release makes it simple to raise/lower cot without assuming full weight of cot and body
- Load end legs unfold and lock automatically when unloading the cot and fold again when loaded into vehicle
- Can act as 2-wheeled stretcher when undercarriage is folded and locked making it highly maneuverable in tight areas
- Slides down stairs with built-in stair guides
- Multi-level feature allows the body to be moved easily to and from cot
- Quick-release restraints hold body in place in all areas
- Waterproof mattress makes it easy to clean
- 5" back end tires allow cot to roll in all areas
- Available in OS 900 Lb. Capacity
- Load Cap. 650 lbs (295 kg)

Available with rails as well as tall models for Suburban type vehicles

Part No	Description	Color	Qty/Cs	Lbs/Cs	Kg/Cs
FCALL-COT-100A	78 x 21" Multi Level FirstCall™ Mortuary Cot	Burg	1	62	28
FCALL-COT-100A-OS	78 x 26" Multi Level FirstCall™ Mort Cot - OS	Burg	1	80	36



FirstCall™ Church Truck

Features 6" wheels and three locking frame positions
 Extends to 48" length
 Constructed of anodized rectangular aluminum tubing
 Max load capacity : 1,000 Lbs (454 kg)
 Available in Silver, Gold or Bronze
 Structural 2 year warranty

Part No	Description	Qty/Cs	Lbs/Cs	Kg/Cs
FCALL-CHTRK-300	4' FirstCall™ Church Truck	1	40	18.2




Medicalproducts Ltd. Eagle Pass TX 78852

Toll Free: (800) 345-2922 Fax: (800) 372-5649 Ph: (678) 413-9200

Internet Site: www.medicalproducts.com Email: sales@medicalproducts.com

A Service Disabled Veteran Owned Company



Aluminum Folding Backboard

Unit features four nylon patient restrainer straps
 Constructed of anodized rectangular aluminum tubing
 Max load capacity : 400 Lbs (181 kg)
 Dimensions when folded: 36" x 16" x 3"
 Structural 2 year warranty

Part No	Description	Qty/Cs	Lbs/Cs	Kg/Cs
FCALL-JSA-362A	16 W X 72" L FOLDING ALUMINUM BACKBOARD- 400 LB CAP – W/O PEG FEET	1	19	8.6
FCALL-JSA-362A-P	16 W X 72" L FOLDING ALUMINUM BACKBOARD- 400 LB CAP – WITH PEG FEET	1	20	9



Aluminum Folding Pole Stretcher

Easy set up and compact for storing
 Medical Corps type light weight litter with wood handles
 Max load capacity : 400 Lbs (181 kg)
 Dimensions when folded: 43 1/2" x 3 1/2" x 7 1/2"
 Cover is durable 18 oz vinyl coated nylon

Part No	Description	Qty/Cs	Lbs/Cs	Kg/Cs
FCALL-JSA-601-NA	22 W X 87" L FOLDING ALUMINUM POLE STRETCHER 400 lb CAP	1	19	8.6



Rugged Plastic Backboard

Rugged one piece lightweight polyethylene spine board
 12 large hand holes and built in runners for easy handling
 Max load capacity : 400 Lbs (181 kg)
 Low profile and X-ray translucent
 Bright yellow color

Part No	Description	Qty/Cs	Lbs/Cs	Kg/Cs
FCALL-JSA-365	16" W X 87" L Plastic Backboard LB 400 lb CAP	1	19	8.6



Medicalproducts Ltd. Eagle Pass TX 78852

Toll Free: (800) 345-2922 Fax: (800) 372-5649 Ph: (678) 413-9200

Internet Site: www.medicalproducts.com Email: sales@medicalproducts.com

A Service Disabled Veteran Owned Company